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```
# Imports
import pylink
import platform
from EyeLinkCoreGraphicsPsychoPy import EyeLinkCoreGraphicsPsychoPy
from string import ascii_letters, digits
import eyetracking

# Set to False if you want to use eyetracking
dummy_mode = True

# Set name for data folder for EDFs
eyetracker_data_folder = eyetracking.create_eyetracking_folder()
# Creates by default a folder with the name "eyetracker data"
```

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```
# Variable for naming EDF
edf_fname = expInfo['participant']
# States that EDF name should be the same as the behavioral experiment

# Check if EDF-name is valid
eyetracking.edf_name_check(edf_fname)
# => should not be longer than 8 characters and only contain letters, digits and
# underscores

# Connects host PC and stimulus PC
el_tracker = eyetracking.connect_tracker(dummy_mode)

edf_file = eyetracking.open_edf_file(el_tracker, edf_fname)

eyetracking.configure_tracker(el_tracker, dummy_mode)
```

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```
eyetracking.calibration_validation(el_tracker, win, dummy_mode)
```

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```
eyetracking.start_recording(el_tracker)
```

```
eyetracking.stop_recording(el_tracker)
```

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```
eyetracking.if_force_end_while_recording(el_tracker, edf_file, edf_fname,  
eyetracker_data_folder)
```

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